



The Alpha Turning Disc can make transfers between bed and (wheel)chair easier for clients that are able to stand. The start and stop function of the Alpha Turning Disc allows movements to be precisely controlled. We recommend the use of a Turning Disc in combination with the Alpha Support Belts.

Technical specification

Product Name	Material	Colour	Max Lifting weigth	Guarantee	Storage/transport
Turning Disc	ABS coated plastic fabric, latex free	beige	200 kg	3 Jahre	5 °C – 40 °C protect against sunlight

Item no.	Dimensions in cm	Weight	UDI
PM-5080	Ø 40	1350 g	(01)04250355350800(10)LOT

Hygiene and maintenance

The product can be wet cleaned. Use warm water and a mild cleaning agent. If applicable, wipe the product off with a moist cloth and then let it dry thoroughly. Don't use any scouring or bleaching agents or any sharp-edged objects for cleaning. When necessary, we recommend wipe disinfection with the product Trionic® or comparable products. The use of bleaching agents or alcohol can permanently impair the properties of the product.



Made in Germany